

## Technical Data Sheet

Article-No.: **09271-M**  
Brand / **TOP-COMFORT**  
Description: **Urinal for men**



SIZES		
length	width	height
310 mm	85 mm	0 mm

PRODUCT DESCRIPTION							
material	<input type="checkbox"/> Polypropylene (PP)	<input checked="" type="checkbox"/> Polyethylene (PE)	<input type="checkbox"/> PP non woven	<input type="checkbox"/> PP non woven with PE coating	<input type="checkbox"/> paper	<input type="checkbox"/> pulp	PVC
colour	<input type="checkbox"/> blue [B]	<input type="checkbox"/> green [G]	<input type="checkbox"/> purple [L]	<input type="checkbox"/> lightblue [LB]	<input type="checkbox"/> lime [Lime]	<input type="checkbox"/> orange [O]	<input type="checkbox"/> black [SCH]
	<input type="checkbox"/> red [R]	<input type="checkbox"/> rose [ROSE]	<input checked="" type="checkbox"/> white [W]	<input type="checkbox"/> yellow [Y]	<input type="checkbox"/> silver	<input type="checkbox"/> natur	<input checked="" type="checkbox"/> transparent
characteristics	<input type="checkbox"/> single use	<input type="checkbox"/> double bladed	<input type="checkbox"/> 1 ply	<input type="checkbox"/> lint free	<input checked="" type="checkbox"/> with white lid	<input checked="" type="checkbox"/> Filling quantity 1.000 ml	
description	Disinfectable and autoclavable up to 60°C						

REGULATORISCHE ANFORDERUNGEN							
PPE-Regulation (EU) 2016/425	<input type="checkbox"/> Category I	<input type="checkbox"/> Category II	<input type="checkbox"/> Category III	<input checked="" type="checkbox"/> no PPE-article			
MD-Regulation (EU) 2017/745	<input type="checkbox"/> Class I	<input type="checkbox"/> Class II	<input type="checkbox"/> Class III	<input type="checkbox"/> sterile	<input type="checkbox"/> measuring function	<input checked="" type="checkbox"/> no medical device	
Food Contact (EG) 1935/2004	<input type="checkbox"/> acidic foods	<input type="checkbox"/> aqueous foods	<input type="checkbox"/> fatty foods	<input type="checkbox"/> alcoholic foods	<input type="checkbox"/> dry foods	<input checked="" type="checkbox"/> not approved for food-contact	

## Technical Data Sheet

Article-No.: **09271-M**  
Brand / **TOP-COMFORT**  
Description: **Urinal for men**



LOGISTIC DATA SUBPACKING		
general information		
material	bag	
pieces per subpacking	1	
EAN subpacking blue [B]	-	
EAN subpacking green [G]	-	
EAN subpacking purple [L]	-	
EAN subpacking lightblue [LB]	-	
EAN subpacking lime [LIME]	-	
EAN subpacking orange [O]	-	
EAN subpacking black [SCH]	-	
EAN subpacking rot [R]	-	
EAN subpacking rose [ROSE]	-	
EAN subpacking white [W]	-	
EAN subpacking yellow [Y]	-	
EAN subpacking silver	-	
EAN subpacking natur	-	
EAN subpacking	-	
measures & size		
length	470 mm	
width	180 mm	
height	0 mm	
weights		
size	net weight	brut weight
uni	100 g	105 g

LOGISTIC DATA OUTER PACKING		
general information		
material	carton	
subpackings per outer packing	25	
EAN outer packing blue [B]	-	
EAN outer packing green [G]	-	
EAN outer packing purple [L]	-	
EAN outer packing lightblue [LB]	-	
EAN outer packing lime [LIME]	-	
EAN outer packing orange [O]	-	
EAN outer packing black [SCH]	-	
EAN outer packing rot [R]	-	
EAN outer packing rose [ROSE]	-	
EAN outer packing white [W]	0	
EAN outer packing yellow [Y]	-	
EAN outer packing silver	-	
EAN outer packing natur	-	
EAN outer packing	-	
measures & size		
length	500 mm	
width	320 mm	
height	420 mm	
weights		
size	net weight	brut weight
uni	2.625 g	3.125 g

LOGISTIC DATA palette		
general information		
kind of palette	Europalette	
measures & size		
cartons per layer	4	
layer per palette	6	
height of the palette	267 cm	
weights		
size	net weight	brut weight
uni	75 kg	100 kg