

FOLDABLE PATIENT STRETCHER



DESCRIPTION

The foldable patient stretcher is suitable for patient transport and consists of a light aluminium frame and a robust PVC lying surface with two safety belts. Thanks to the carrying straps, the patient can be reliably secured. For storage, the patient stretcher can be folded in the middle thanks to the folding mechanism and stored in the carrying bag provided to save space. When folded, the stretcher takes up very little space and can be easily stored in first aid cabinets or vehicles.

PRODUCT DETAILS

- Foldable patient stretcher
- With PVC lying surface and aluminium frame
- With wrapped carrying handles
- Two safety straps ensure safe transport
- Can be folded up
- Stretcher with 4 small feet
- Can be stowed away to save space
- Load capacity up to 150 kg
- Dimensions when unfolded: 208 x 55 x 13 cm
- Dimensions when folded: 104 x 17 x 9 cm
- Carrying bag included
- 1 piece

Overview




Article Number:
T135444

Category:
Emergency Equipment

EAN:
4260547365819

Teqler is a trademark of NetMed S.à.r.l.

 NetMed S.à.r.l.
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  info@netmed.lu

Social Media

