



# Typhoid IgG/IgM Rapid Test Kit (Colloidal Gold Method)

The Typhoid IgG/IgM Rapid Test Kit is a lateral flow immunoassay assay designed for qualitative detection and differentiation of IgG and IgM Typhoid antibodies produced against Salmonella Typhi (*S. typhi*) in human whole blood, serum or plasma. This easy-to-use format with rapid results allows point-of-care testing, providing information to assist with treatment and patient management decisions.

## Test Procedure

- **Specimen:** Whole blood, serum or plasma.
- **Specimen collection:** Specimens of all types are collected by conventional methods.
- **Test method:**



## Specifications

Components	Test card/Sample diluent/ Sample Droppers
Pack Size	25Tests/Kit; 50 Tests/Kit
Storage Conditions	2°C~30°C ;Valid for 24 months
Test samples	Whole blood, Serum or Plasma
Reaction conditions	15~30°C 15~20 minutes
Special instrument	No need

## Product Informations

Product Name	Cat.No	Package	Sample Diluent	Sample Droppers	Sample Type
Typhoid IgG/IgM Rapid Test Kit (Colloidal Gold Method)	BSK11S1S	25 Tests/Kit	1 × 3 mL/Tube	25	Whole blood, serum & plasma
	BSK11M1S	50 Tests/Kit	2 × 3 mL/Tube	50	